

# KERN 313-066

**KERN**

Weight set, E2, 1 mg - 1 kg, Knob, stainless steel polished, in aluminium case



Gross weight approx. 3,6 kg

Shipping weight 3,5 kg

## Optional cases

Case - Aluminium article number 313-060-600

## Pictograms

## STANDARD



## Category

|                  |            |
|------------------|------------|
| Brand            | KERN       |
| Product category | Weight set |
| Product group    | OIML E2    |
| Product family   | 313        |

## Measuring System

|                                |             |
|--------------------------------|-------------|
| Conformity assessment possible | ✓           |
| Test weight - nominal value    | 1 mg - 1 kg |
| OIML-compliant                 | ✓           |
| Weight OIML class              | E2          |

## Construction

|                    |                          |
|--------------------|--------------------------|
| Dimension (W×D×H)  | 230×140×110 mm           |
| Shape weight       | Knob                     |
| Material           | stainless steel polished |
| Case of weight set | in aluminium case        |

## Services (optional)

|                                              |         |
|----------------------------------------------|---------|
| Article number for conformity assessment     | 952-305 |
| Article number for DAkkS calibration         | 962-305 |
| Article number for guaranteed plus-deviation | 921-305 |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Delivery time                | 1 d            |
| Dimensions packaging (W×D×H) | 230×140×110 mm |
| Net weight                   | 2,7 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 2,8 kg         |