

KERN 326-36

KERN

Individual weight, OIML Class F1, 5 kg, Block, stainless steel polished



Dimensions packaging (W×D×H)	150×68×96 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

Optional cases

Case - Plastic article number	nicht verfügbar
Case - Aluminium article number	346-060-600
Case - Wood article number	nicht verfügbar

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Conformity assessment possible	✓
Test weight - volume uncertainty (cm³)	7,8125
Test weight - Density (kg/m³)	7950
Test weight - Density uncertainty (kg/m³)	140
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Test weight - nominal value	5 kg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	25 mg

Construction

Shape weight	Block
Dimensions weight (W×D×H)	150×68×96 mm
Material	stainless steel polished

Services (optional)

Article number for DAkkS calibration	962-443
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
---------------	-----