

KERN 326-39

KERN

Individual weight, OIML Class F1, 50 kg, Block, stainless steel polished



Dimensions packaging (W×D×H)	295×158×191 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Optional cases

Case - Plastic article number	nicht verfügbar
Case - Aluminium article number	346-090-600
Case - Wood article number	nicht verfügbar

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

Measuring System

Conformity assessment possible	✓
Test weight - volume uncertainty (cm ³)	78,125
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Test weight - nominal value	50 kg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	250 mg

Construction

Shape weight	Block
Dimensions weight (W×D×H)	295×158×191 mm
Material	stainless steel polished

Services (optional)

Article number for DAkkS calibration	962-446
--------------------------------------	---------

Packing & Shipping

Delivery time	2 d
---------------	-----